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规格: 100mL/500mL

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**Cat:** G3005

**Size:** 100mL/500mL

**Storage:** 2-8°C, avoid light, valid for 1 year.

### Introduction

2,3,5-triphenyltetrazolium chloride (TTC) is a lipid soluble light sensitive complex, which is commonly used to detect the viability of seeds and the ischemic infarction of mammalian tissues. It is the proton receptor of pyridyl-nucleoside structural enzyme system in respiratory chain, which can react with dehydrogenase in normal tissues and show red color. The dehydrogenase in the living cells of the viable seed embryo tissue can reduce TTC to insoluble red stable triphenylmethane (TTF). If the embryo dies or the viability of the embryo declines, it cannot be stained or stained shallowly. Therefore, the viability of seeds can be determined according to the location and depth of staining.

### Protocols(*for reference only*)

1. According to the experimental design, make the tissue model and slice.
2. Dyeing: soak the tissue or slice in TTC Solution, 2% and put it into a constant temperature box at 30-37 °C for 15min-30min.
3. After dyeing, wash with water and fix in formalin or polyformaldehyde.
4. Observe the staining result and photograph.

### Result

Normal tissues	Red
Cerebral infarction area	Not colored

### Note

1. TTC Solution, 2% can not be reused if it turns red.
2. The section should be absolutely immersed in TTC Solution, 2% when staining.
3. Immediately observe after staining with TTC Solution, 2% and the section can not be placed for a long time.
4. For your safety and health, please wear experimental clothes and disposable gloves.